



# ETHYLENE OXIDE GAS CARTRIDGES ETİLEN OKSİT GAZ KARTUŞLARI



**% 100 EO**

Axis provides Ethylene Oxide in aluminium cartridges with controlled quantities and packed under strict safety regulation that meet the requirements of the variety of existing equipment. Ethylene Oxide of high pureness. Single-use cartridges. High effectiveness as a bactericide, fungicide and antiviral. Compatible with any sterilizer in the market.

Tek kullanımlık alüminyum gövdeli Etilen Oksit kartuşlarıdır. %100 konsantrasyonda olup sterilizasyon güvencesini maksimuma çıkarmaktadır. Sınıfında her marka cihazda kullanıma uygundur.

CODE	CAPACITY / KAPASİTE		
	Diameter Çap	Height Yükseklik	Net Content Muhteviyat
AL-7	25 mm	85 mm	10 g
AL-25	25 mm	105 mm	30 g
AL-67	35 mm	105 mm	67 g
AL-100	35 mm	165 mm	100 g
AX-67	38 mm	105 mm	67 g
AX-100	38 mm	165 mm	100 g
AX-127	38 mm	165 mm	127 g
AX-134	38 mm	165 mm	134 g
AQ-70	45 mm	178 mm	170 g
AQ-120	50 mm	133 mm	120 g
AQ-130	50 mm	133 mm	130 g
AQ-200	50 mm	174 mm	200 g
AQ-230	50 mm	174 mm	230 g

**@x** **AXIS**

[www.axis.com.tr](http://www.axis.com.tr)